

STAPLER

SKIN STAPLER SYSTEM

The benefits are widely known and accepted in human surgery: better apposition of skin edges with better cosmetic results, less reactivity of stainless steel staples and significantly faster wound closure with staples.

How to use the system

First approximate the skin edges in an everted fashion to desired height.

Gently place the nose of the instrument against the middle of the incision. Squeeze the trigger completely and release (failure to completely squeeze the trigger can result in incomplete staple formation and insufficient wound approximation). Repeat this procedure, taking care to space the staples evenly, until the incision is completely closed. Check to ensure staples are correctly formed (any malformed staples should be removed and replaced with a new staple).

Y4063 SMI SKIN STAPLER

This skin stapler can be used for a variety of skin closures. It is ergonomic and user friendly, with a clear view of operating site and is easy to use in checking remaining staples. Contains 35 staples (6,4 x 4mm)



Y4063

Y4061 SURGICLOSE™

HANDLE/APPLICATOR

The handle is reusable and disposable which allows for reloading of sterile multi-feed staples. Extremely shallow tissue channel and eversion assisting design.



Y4061

Y4062 SURICLOSE STAPLE CARTRIDGES
Staple cartridges are especially designed to reduce the chance of staples jamming. Skin staple after closure is 6.65mm wide x 4.65mm high, a better length for thicker animal skin versus thin human skin. The 20 staples provided per staple cartridge are sufficient for a 100-125mm incision.



Y4062

Y4060 STAPLE REMOVER
This versatile staple remover will remove any brand of surgical staples, simply and efficiently. The staple removal tool removes staples by the tip of the remover under each staple, which saves time and provides reliable, consistent performance.



Y4060

Y4098 SURGIBOND
A N-butyl cyanoacrylate tissue adhesive for:

- ✓ Wound closure
- ✓ For additional suture support
- ✓ Stops bleeding
- ✓ Protects against infections
- ✓ Good visibility due to blue colour
- ✓ Acts within seconds
- ✓ No unnecessary sutures or dressings
- ✓ Indicated for: lacerations, abrasions, suture/ staple line sealing, tooth extractions & oral surgeries



Y4098