


COLLARS

BIRD COLLARS

Transparent disposable collar for birds. It prevents birds from picking at feathers and wounds. Easy to apply with an office stapler.

Sizes (refer diagram below, 1):

- Y2640 X-small (8cm)
- Y2641 Small (10cm)
- Y2642 Medium (14cm)
- Y2643 Large (20cm)



Y2640-Y2643

CLASSIC E - COLLARS

Quick and simple “clac” closing system. Fully transparent material. The newly developed material will allow the animal all-round vision to avoid walking into furniture, etc.

- ✓ Soft and flexible.
- ✓ Easy to apply to patient.
- ✓ Made from environmentally-friendly 100% polypropylene.

Sizes (refer diagram below):

Classic E-Collar		Fitting for neck Ø
Y2625	7,5cm	6-8cm
Y2626	10cm	7-8cm
Y2627	12,5cm	9-11cm
Y2628	15cm	10-11cm
Y2629	20cm	12-14cm
Y2630	25cm	14-16cm
Y2631	30cm	15-18cm
Y2633	35cm	19-22cm
Y2632	40cm	22-25cm



Y2625-Y2632

