

SMALL ANIMAL MOUTH GAGS

GREY MOUTH GAG

Spring action with nylon inserts for reduced trauma.

Sizes:

- Y1928 Feline (7,5cm)
- Y1925 Small (11,0cm)
- Y1926 Medium (15,0cm)
- Y1927 Large (17,0cm)

Spare rubber stopper available on request.



Y1925-Y1928

SMITH BAXTER MOUTH GAG

Sizes:

- Y1934 Feline size
 - 1 = 5cm
 - 2 & 3 = 3cm
 - 4 = 3cm
- Y1930 Small (cats and very small dogs)
 - 1 = 5cm
 - 2 & 3 = 6cm
 - 4 = 3cm
- Y1931 Medium (small and medium dogs)
 - 1 = 6cm
 - 2 & 3 = 9cm
 - 4 = 4cm
- Y1932 Large (medium and large dogs)



Y1930-Y1934