


SURGICAL

Y6111 BONE MARROW BIOPSY NEEDLE (16G) [JORVET]

Ergonomically designed, with a large handgrip for increased torque. Adjustable depth stop from 10mm - 37mm shaft. Luer lock syringe attachment. Similar to Jamshidi style. Can also be used for intraosseus fluid administration. *For tips on how to use this device, please refer to Section Y — see page 136.*




Y6111

DISPOSABLE SKIN BIOPSY PUNCH

Sterile biopsy punch for sampling of skin tissue.

Sizes:

- Y2575 4mm
- Y2574 6mm



Y2575 • Y2574


TRU-CORE II AUTO BIOPSY NEEDLE (Not Autoclavable)

The next generation of soft tissue biopsy instruments (sterile). The Tru-Core II is widely used in human medicine because:

- ✓ Ease-of-use
- ✓ Better, more representative samples
- ✓ Uses 1, 2, 3 process clearly labeled on the instrument
- ✓ Echogenic tip for use with an ultrasound

Sizes:

- ✓ 14G x 10cm
- ✓ 14G x 16cm
- ✓ 16G x 10cm
- ✓ 16G x 16cm
- ✓ 18G x 10cm



TRU-CORE II AUTO BIOPSY NEEDLE