

## SECTION M – MOUTH GAGS: LARGE ANIMALS

### BOVINE

#### Y1910 BAYER MOUTH GAG

To be used in cattle when passing a stomach tube or examining the oral cavity.

#### DRINKWATER MOUTH GAG

An efficient and reliable gag for use in cattle when passing a stomach tube or examining the mouth. When properly fixed between the upper and lower molar teeth, is self-retaining.

#### Y1901 Left gag

#### Y1900 Right gag



Y1910



Y1900 • Y1901



Y1915

#### Y1915 GUNTHER HEPKE MOUTH GAG

For use in equines when passing endotracheal tubes.

#### Y1917 CLIMAX MOUTH GAG


Forged stainless steel complete with halter — ideal for dental procedures.



Y1917

#### Y1921 HHP-MODEL EQUINE MOUTH GAG

Made of solid stainless steel with a moulded rubber block. The rubber block is made of very durable material and is atraumatic so that tooth damage does not occur. In addition the wedge shape and the compliance of the rubber cause less stress to the horse. The rubber block is exchangeable — when it wears out it can easily be replaced.



Y1921

#### Y1920 SCHOUPE MOUTH GAG

Reduces the danger of slipping since it is held in place by the equine cheek area. Inserted without difficulty. It permits greater freedom in treating any part of the oral cavity.



Y1920