

BOVINE

ARMLENGTH (RECTAL) GLOVES (NON-STERILE 100'S)



Y2001

Y2001 SMI green soft gloves are thin yet strong, disposable gloves.

Y2002 Clear, super sensitive disposable gloves.

New Green shoulder glove



Y2002



Y2802

Y2802 ARTIFICIAL VAGINA (AV)
Complete AV for cattle. Consists of a rubber AV, latex liner and cone.
Available on request.

Y2800 AV CONES

Size: 19cm



Y2800



Y2801

Y2801 AV LINERS

Size: 60cm

Y2805 BLUNT EYEHOOK

Size: 6cm



Y2805



Y2811

CALF PULLER

This instrument is used during the calving process when the cow is unable to deliver the calf by herself. The ropes are attached to the calf's legs and the calf is extracted using the calf puller.

Y2811 Calf puller ropes (a pair of replacement ropes).

Y2810 Calf puller with set of angled pieces to keep from slipping.

Y2809 Calf puller with support holding bows hooking over/around hips to prevent slipping.



Y2810



Y2809

CALVING CHAIN

Stainless steel calving chain.

Length:

Y2814 80cm with 1 flat inset link

Y2815 190cm with 2 flat inset links



Y2815

SECTION N – OBSTETRICAL EQUIPMENT

Y2816 CALVING CHAIN HANDLE

Y2820 CALVING ROPE SET (HAUPTNER)
Triple set calving ropes with head rope (99cm) and 2 leg ropes (115cm each).

Y3550 COLOSTRUM DENSIMETER
For quick evaluation of IgG content in bovine colostrum. The Densimeter measures the specific gravity (SG) of the colostrum. Colostrum with a low SG has a low content of IgG.

ELECTROJAC 6

Electrojac 6 is a fully portable ejaculation unit that will operate manually or automatically through 32 increasing power output cycles. The dwell switch eases semen collection, especially from older bulls that may be difficult to stimulate.

Features:

- ✓ Improved output curve for more efficiency when collecting
- ✓ Early cycles more comfortable for younger bulls and later cycles have 5% more power for older bulls
- ✓ Long-lasting, powerful, gel cell battery (The unit can be left plugged-in indefinitely after each use to ensure a fully charged battery whenever the unit is needed).
- ✓ Unparalleled simplicity of use
- ✓ Unique dwell switch for cycle pausing (To control automatic cycles: when the dwell switch is depressed & held, the Electrojac 6 will remain at that cycle [2 seconds on, 2 seconds off] until the switch is released.)
- ✓ Bright numeric display (In automatic mode, the numerals indicate which of the 32 cycles the unit is in. In manual mode, the numerals will blink and indicate output voltage.)



Y2816



Y2820




Y3550



Y2835

ElectroJac 6 consists of:

- ✓ Electronics instrument case
- ✓ Dwell cable assembly
- ✓ Probe cable assembly
- ✓ Semen collection director cone & vial
- ✓ Probe of choice
- ✓ Carrying case
- ✓ AC battery charger
- ✓ Metric scrotal tape
- ✓ Transformer (provided for SA)
- ✓ Remote Control



Y2835/I

ELECTROJAC 6

Y2835 with 76mm (3") bull probe

Y2835/I with 64mm (2,5") bull probe

Extra probes available on request
64mm / 51mm Bull probe
25mm Ram probe

Y2835/I METRIC SCROTAL TAPE
Available separately from the complete set.



Y2831



Y2830

EMBRYOTOMY

Y2831 EMBRYOTOME HANDLES (pair)

EMBRYOTOMY WIRE

Used for embryotomy and dehorning.

Y2830 Jorvet embryotomy wire (10m)

Y2829 Hauptner embryotomy wire (10,8m)

Y2826 KREY SCHÖTTLER HOOK

This double action hook will bite into part of the foetus in the process of an embryotomy. Should the instrument slip, the hooks will close automatically and prevent injury.

Length: 20cm

Y2825 LINDE'S EMBRYOTOMY KNIFE

One large finger ring and handle adapted to the palm.



Y2829



Y2826



Y2825

SECTION N – OBSTETRICAL EQUIPMENT

THYGESEN EMBRYOTOME

Two-tube model with an adjustable handle for calving chain. A complete embryotome consists of:

- Y2840 A nickel plated embryotome (100cm)
- Y2842 A one-piece threader (110cm)



Y2840 • Y2842

- Y2861 NITROGEN FLASK FORCEPS
Long narrow forcep to take AI straws out of the nitrogen flask. 25cm



Y28561



Y2822

OBSTETRIC WIRE GUIDE

This instrument will bring wire around any tight or difficult location on the foetus.

Types:


- Y2821 Light wire guide
- Y2822 Heavy wire guide

- Y3551 UMBILICAL TAPE
Braided, non-sterile polyester tape in plastic bottles.
Size: 6mm x 18m



Y3551

- Y3520 UTERINE INFUSION PIPETTES
High quality non-sterile PVC pipettes with rounded polished insertion end to prevent tissue trauma. Large bore with flexible syringe adapter for infusing viscous fluids. Fits onto catheter tip of dosing/flushing syringe (Y4415). Packaged 25 per poly bag.
Length: 45cm



Y3520

VAGINAL SPECULA

STAINLESS STEEL WITH TWO BLADES (“duckbill type”)

Sizes:

- Y3515 27cm
- Y3516 35cm



Y3515 • Y3516




Plexiglass Speculum

PLEXIGLASS SPECULUM WITH LIGHT HANDLE (Hauptner)

Size: 350 x 40mm. Available on request.

VULVA PROLAPSE SET

Y3530 For repair of bovine uterine or vaginal prolapse. A set consists of two plastic spacers and three prolapse pins with wooden balls.



Y3530

VULVA SUTURE NEEDLE

Ideal for use with vulva sutures. Can also be used for piercing the skin for rumenotomy sutures.

Y2845 Bühner needle (30cm)

Y2846 Gerlach needle (19cm)



Y2845



Y2846

CANINE

Y2860 DANISH WHELPING FORCEPS

VAGINAL SPECULA

CUSCO STAINLESS STEEL WITH TWO BLADES (“duckbill type”).

Available on request.

Y2865 PLEXIGLASS SPECULUM WITH LIGHT HANDLE (Hauptner)

Size: 120 x 20mm

Y2865/T Perspex tube only



Y2860



Cusco Speculum

EQUINE

AV SET (MISSOURI TYPE)

Model Missouri, complete with liner, leather cover and two collecting bottles. The vagina is to be filled with water.

Length: 46cm.

Available on request.




Y2865



AV Set

MARE UTERINE CATHETERS

- Y1270 EQUIVET uterine flushing catheter:
Sterile 190cm catheter with bore spikes. Ensures easy and secure coupling of flushing tube and flushing liquid.



Y1270

VAGINAL SPECULA

- Y3510 POLANSKY SPECULUM
Stainless steel with 3 blades.
Length: 35cm.
Available on request.



Y3510



Y3500

- Y3500 PLEXIGLASS SPECULUM WITH LIGHT HANDLE (Hauptner)
Size: 420 x 50mm

OVINE

UTERINE SUPPORT EWE

Plastic support to prevent prolapse around the time of lambing. Once inserted into the vagina it can be tied to the fleece. Sterilise and use a lubricant before inserting.
Pack of 3.



Uterine Support Ewe

VAGINAL SPECULA

- Y3518 STAINLESS STEEL WITH TWO BLADES (“duckbill type”)



Y3518



Y3501

- Y3501 PLEXIGLASS SPECULUM WITH LIGHT HANDLE (Hauptner)
Size: 250 x 30mm