

CATTLE

Y3315 ANTI-KICK BAR
For immobilising cattle.

BULL NOSE HOLDER

Y3311 HARMS MODEL
With slide adjuster, stainless steel.
Length: 19cm

Y3310 TONG TYPE

Y3320 PRODDER (KAWE)
Sturdy prodder with a strong shock.
Length: 22,5cm (Extension handle
available on request).



Y3315



Y3311



Y3310



Y3320

EQUINE

Y3312 NOSE TWITCH (EQUIVET)
Aluminium with one metre nylon rope and snap hook.

SMALL ANIMAL

DOG MUZZLES

Made of nylon, it fits snugly over the dog's muzzle and fastens around the neck with a quick release clip.

Sizes:

Y2660 Small

Y2661 Medium

Y2662 Large

Y2663 X-large

GRABBERS

Y 2652 CAT GRABBER
This lightweight grasper was designed primarily for the effective seizure and control of cats. It is ideal for reaching into the back of a cage to secure a difficult cat. The smoothly curved and finished jaws are plastic-covered and designed to minimize discomfort. If necessary the animal can be quickly released.
Length: 76cm



Y3312



Y2660-Y2663



Y2652


DOG GRABBER

Manufactured from aluminium. Strong plastic covered steel cable to capture the animal, can be locked in position by means of a twist action of the handle. Spring mechanism allows instant release.

Lengths:

Y2651 100cm

Y2650 150cm



Y2650 • Y2651

SAFETY STICK POLE SYRINGE

This syringe allows you to obtain a proper needle stick, then push the dosage plunger on the operator end. It is constructed of aircraft aluminium alloy stainless steel. Black anodised satin finish.

Features:


- ✓ Uses a standard Monoject 3cc syringe.
- ✓ Mechanism cannot be damaged by even the most powerful animal bite.
- ✓ Easy to inject through cage or from a safe distance.

Lengths:

Y2653 914mm

Y2654 1,2m

Y2655 1,8m




Y2653 • Y2654 • Y2655

Y2670 SOFT CARRY STRETCHER

Great for transporting injured animals. Absence of rigid poles allows carrying patients in tight corners. Six padded handles for carrying comfort.

Washable nylon.

Size: 75 x 125cm



Y2670