

Y1860 COTTON WOOL HOLDER

## DRESSING TROLLEYS (with rail)

*Available on request.*

## DRUG CABINETS

Small (Lockable) – 37 x 24 x 46cm

Large (Lockable) – 31 x 46 x 61cm

Deadly dangerous cabinet.

*Available on request.*

## INSTRUMENT TRAY AND LID

Sizes:

Y1830 200 x 125 x 50mm

Y1820 250 x 140 x 50mm

Y1821 325 x 265 x 50mm

Y1822 300 x 190 x 50mm

Y1823 420 x 300 x 50mm

## INSTRUMENT TROLLEYS

(without rail)

Stainless steel shelves and epoxy

coated frame. Height: 85cm

Sizes:

Y1810 450 x 450mm

Y1811 610 x 450mm

Y1812 910 x 450mm

## KICK-ABOUT with:

Y1835 Bucket

Y1836 Bowl



Y1860



Drug Cabinet



Y1830 • Y1820-Y1823



Y1810-Y1812



Y1835



Y1836

## KIDNEY DISHES (stainless steel)

Sizes:

- Y1800 15cm
- Y1801 20cm
- Y1802 25cm
- Y1803 30cm



Y1800-Y1803

## Y1851 MAYO TABLE

Stainless steel tray and epoxy coated frame

Tray size: 60 cm (L) x 40 cm (W)

Height: 96cm



Y1851

## Y1840 MOBILE DRIPSTAND




Y1840

## Y1855 SPAY CRADLE

2-Bar spay cradle for positioning animals for an operation.

(Chrome plated and aluminium).

(A complete surgical instrument spay set is also available – refer page 92).



Y1855