

**ANAESTHESIA**

**ENDOTRACHEAL TUBE CONNECTORS**

Individually packaged Jorvet non-sterile 15mm tracheal tube connectors.

Sizes:

- 2,5mm
- 3mm
- 3,5mm
- 4mm
- 4,5mm
- 5mm
- 5,5mm
- 6mm
- 6,5mm
- 7mm
- 8,5mm
- 9mm
- 9,5mm
- 10mm
- 10,5mm
- 11mm



Anaesthetic Connectors

**ENDOTRACHEAL TUBES**

*SMALL ANIMAL*

**PVC, CUFFED ENDOTRACHEAL TUBES with CONNECTORS**

Blue tinted PVC tubes soften at body temperature, reducing post-operative irritation.

Sizes:

- Y1033 3mm
- Y1034 3,5mm
- Y1035 4,0mm
- Y1036 4,5mm
- Y1037 5,0mm
- Y1038 5,5mm
- Y1039 6,0mm
- Y1040 6,5mm
- Y1041 7,0mm
- Y1042 7,5mm
- Y1043 8,0mm
- Y1044 8,5mm
- Y1045 9,0mm
- Y1046 9,5mm
- Y1047 10,0mm
- Y1048 10,5mm
- Y1049 11,0mm



Y1033-Y1049

## SECTION X — THEATRE EQUIPMENT

SILICONE, CUFFED  
ENDOTRACHEAL TUBES with  
CONNECTORS

CANINE

Sizes:

- Y1065 12mm x 44cm
- Y1066 14mm x 48cm




Y1065 • Y1066

RED RUBBER, CUFFED  
ENDOTRACHEAL TUBES

The standard oral Magill ET tube to which is vulcanised an inflatable cuff and pilot balloon. (ET tube connector/adaptor not included).

Sizes:

- Y1030 2,5mm
- Y1001 3,0mm
- Y1002 3,5mm
- Y1003 4,0mm
- Y1004 4,5mm
- Y1005 5,0mm
- Y1006 5,5mm
- Y1007 6,0mm
- Y1008 6,5mm
- Y1009 7,0mm
- Y1010 7,5mm
- Y1011 8,0mm
- Y1012 8,5mm
- Y1013 9,0mm
- Y1014 9,5mm
- Y1015 10,0mm
- Y1016 10,5mm
- Y1017 11,0mm



Y1001-Y1020 • Y1030

Special sizes with adaptor are available  
on request:

- Y1018 12mm x 48cm
- Y1019 14mm x 48cm
- Y1020 16mm x 48cm

## EQUINE with ADAPTORS/CONNECTORS

(Available on request)

Sizes:


- Y1025 20mm x 112cm
- Y1026 25mm x 117cm
- Y1027 30mm x 117cm
- Y1028 35mm x 117cm
- Y1029 40mm x 122cm

## ENDOTRACHEAL TUBE CLEANING BRUSHES

Nylon bristles attached to a stainless steel shaft. The rigid structure allows for easy cleaning of ET tubes.

Sizes:

- Y1370 Small (5mm x 25cm)
- Y1371 Medium (8mm x 32cm)
- Y1372 Large (10mm x 42cm)



Y1370-Y1372

## LARYNGOSCOPE

### MACINTOSH BLADE (curved)

Sizes:

- Y1085 #0 (80mm)
- Y1086 #1 (92mm)
- Y1087 #2 (100mm)
- Y1088 #3 (130mm)
- #4 (155mm)
- #5 (175mm)



Y1085-Y1088

### MILLER BLADES (straight)

Sizes:

- # 00 (65mm)
- Y1081 # 0 (77mm)
- Y1082 # 1 (102mm)
- Y1083 # 2 (153mm)
- Y1084 # 3 (195mm)
- # 4 (205mm)



Y1081-Y1084

- Y1080 LARYNGOSCOPE HANDLE  
A standard, medium size battery handle.



Y1080

## PLEXIMASK

Anaesthesia mask with low dead space. Clear plexiglas mask with black rubber diaphragm.

Sizes:

- Y1074 Canine (13cm opening)
- Y1073 Feline (9cm opening)
- Y1077 X-small (3cm opening)



Y1073-Y1074

## PLEXIMASK REPLACEMENT DIAPHRAGMS

Black rubber diaphragm to fit Plexiglas mask.

- Y1071 Canine
- Y1072 Feline



Y1071



Y1072

## REBREATHING BAGS

Black anti-static rubber rebreathing bag with a standard size ferrule neck. It may be cleaned with a mild detergent and dried at room temperature.

Sizes:

- Y1050 0,5L
- Y1052 1L
- Y1054 2L
- Y1055 3L
- Y1056 4L
- Y1057 5L
- Y1058 6L
- Y1051 15L
- Y1059 30L



Y1050-Y1059

**JORVET RESUSCITATION BAGS**

These durable PVC ventilation bags can be used to ventilate the animal either through face mask or, preferably, endotracheal tube (2 non-return valves allow you to ventilate with oxygen or air). Extremely handy to have in case of emergencies, when assisted ventilation of the animal is needed (a rebreathing bag can be fitted for asthma drugs in cats).

Sizes:

Y1075 Small (cats & small dogs)

Y1076 Large (large breed dogs)



Y1075 • Y1076

**ANAESTHETIC REBREATHING CIRCUITS**

Sizes:

Y1061 Anaesthetic circuit paediatric without bag (disposable)

Y1062 Anaesthetic circuit adult without bag (disposable)

Y1063 Anaesthetic circuit non-breathing with bag (disposable)



Y1061

Y1060 AYRES T-PIECE NON-REBREATHING SYSTEM WITH SCAVENGER ADAPTER (COMPLETE SYSTEM)

A valveless, semi-closed, non-rebreathing anaesthetic system designed especially for the small animal patient (cats and puppies with less than 6.8kg body weight). It provides minimal dead space and less resistance to flow as compared to a full-circle system. An anaesthetic mask can easily be attached to the metal collar (connects in seconds to any Matrix™ anaesthesia machine).



Y1063



Y1060