

Y2565 MINI-CAUTERISER

A small portable cordless cauteriser excellent for small bleeders. A handpiece containing 2 AA batteries plus four replaceable cautery tips in a handy foam padded carrying case. Delivers a temperature up to 1232°C. Ideal for fieldwork.

Cautery Tips:

- ✓ 1,3cm fine electrode (useful for haemostasis control)
- ✓ 1,3cm loop electrode for wide use (such as small skin tumours)
- ✓ 5cm fine electrode for difficult areas
- ✓ 5cm loop electrode



Y2565

SPECIALISED

ANAL SAC EXCISION KIT

The Anal Sac Excision Kit utilises a gel liquified by heat and injected into the anal sac. The gel is in a 2ml cartridge that is first boiled to a liquid, loaded in a dental syringe and with a special blunt needle injected into each anal sac. The gel hardens and gives the sac a readily defined and easily dissected shape for quick removal.

The kit consists of:

- Y6000 8 x 2ml Anal Sac Gel vials
- Y6001 8 x 19g blunt Anal Sac Needles
- Y6002 Metal breech loading dental syringe



Y6000



Y6001



Y6002

Y6040 OPTIVISOR HEADPIECE The Optivisor is a head-worn magnifying glass. The optical-glass quality lenses are precision ground and polished, and the body of the instrument is made of lightweight, durable, impact-resistant material. The headpiece adjusts to all head sizes. The equally-ground prismatic lenses assist the eyes to converge at a given focal distance. This focal distance depends on the power of the lens being used with the instrument.

Note: Lenses must be ordered separately to the headpiece.

Lens Plates:

	No.	Magnifies	Focal distance
Y6041	2	1,5x	50cm
Y6042	3	1,75x	35cm
Y6043	4	2x	25cm
Y6044	5	2,5x	20cm
Y6045	7	2,75x	15cm
Y6046	10	3,5x	10cm

Y6047 VISORLIGHT
With battery pack & cord. Fits onto the OPTIVISOR (sold separately) as an extra light source.



Y6040



Y6047