

### Calibration Table for Shiotz Tonometry

Recommended weight (in both cats and dogs): 7.5g

VALUES IN CATS				VALUES IN DOGS		
Schiotz Scale Reading	IOP (mm Hg) 5.5g wt.	IOP (mm Hg) 7.5g wt.	IOP (mm Hg) 10.0g wt.	IOP (mm Hg) 5.5g wt.	IOP (mm Hg) 7.5g wt.	IOP (mm Hg) 10.0g wt.
0.5	44	73		46	61	75
1.0	42	71		44	59	73
1.5	40	68		43	56	70
2.0	37	65	80	40	53	66
2.5	33	61	76	33	47	61
3.0	30	56	71	26	40	55
3.5	27	48	66	23	35	49
4.0	25	42	61	21	32	44
4.5	24	37	56	20	29	41
5.0	22	34	51	19	27	38
5.5	21	31	47	18	26	36
6.0	20	29	44	17	24	33
6.5	18	27	40	16	23	31
7.0		25	37	15	22	30
7.5	17	24	35		20	28
8.0	16	22	33	14	19	27
8.5	15	21	31	13		25
9.0	14	20	29		18	24
9.5	13	19	27	12	17	23
10.0		18	25		16	22
10.5		17	23	11	15	21
11.0	12	16	22			20
11.5	11	15	20	10	14	19
12.0		14	19		13	18
12.5	10	13	18			17
13.0		12	17		12	16
13.5	9		15	8	11	15
14.0		11	14			
14.5	8	10	13		10	14
15.0			12	7		13
15.5		9	11		9	12
16.0	7	8	10			
16.5			9	6	8	11
17.0	6	7	8			10
17.5		6	7		7	
18.0			6	5		9
18.5	5	5	5		6	
19.0						8
20.0					5	