

PREVENTATIVE MAINTENANCE PROGRAMME — CARDELL®

This program is designed for monitors to go through a like new process, aiding in maximising their life span. By sending your monitor in for a PM we are able to:

- ✓ Perform internal design revisions and update software
- ✓ Replace the NIBP receptacle along with any other connections that may have become loose preventing air leaks at hose connection
- ✓ Battery replacement for longer use between chargings
- ✓ Test all accessories sent in with the monitor
- ✓ Replacement of temperature connection to prevent intermittent temperature readings
- ✓ Front panel keyswitch and housing replacement as needed returning unit to like new condition
- ✓ NIBP Pneumatic check
- ✓ Monitor Safety Leakage Check/ Hi Pot Test
- ✓ Updated labeling
- ✓ Functional check, QA Testing and calibration of the unit for more accurate readings

If you are interested in our Preventative Maintenance Program, please contact us for current pricing.